

### BD FACSMelody Sorter (5532)

<b>Laser</b>	<b>Detector</b>	<b>Mirror</b>	<b>Filter</b>	<b>*Suggested Fluors</b>
<b>Violet (405nm)</b>	V450		448/45	Pacific Blue, BV421
	V510	500LP	528/45	BV510
	V786		755LP	BV786
<b>Blue (488nm)</b>	FITC	507LP	527/32	FITC, GFP, Sytox Green
	PE	560LP	586/42	PE, TdTomato
	PerCP-Cy5.5	665LP	700/54	PerCP-Cy5.5, 7AAD
	PE-Cy7	752LP	783/56	PE-Cy7
<b>Red (640nm)</b>	APC	660/10	660/10	APC, Alexa647, Sytox Red
	APC-Cy7	752LP	783/56	APC-Cy7

\*List of suggested fluorochromes is not exhaustive